



Applied HPC Seminar Series: Introduction to Containerization for HPC

What is Containerization?

- Think of it as a “portable research environment”
- Contains everything your application needs to run:
 - Operating system files
 - Software dependencies
 - Libraries
 - Configuration files
- Like a lightweight, portable computer within your computer

Why Use Containers in Medical Imaging Research?

- **Reproducibility:** Same environment = same results
- **Portability:** Run anywhere (local machine, HPC cluster, cloud)
- **Version Control:** Track changes in your research environment
- **Collaboration:** Share exact analysis environment with colleagues
- **Dependency Management:** No more “works on my machine” problems

Real-World Example: Medical Image Processing

Before Containers:

```
❌ Install specific CUDA version  
❌ Install specific Python version  
❌ Install TensorFlow with GPU support  
❌ Install specific NVIDIA drivers  
❌ Install medical imaging libraries  
❌ Hope everything works together
```

With Containers:

```
✅ Pull container image  
✅ Run your analysis
```

Popular Container Technologies in HPC

Singularity/Apptainer

- Designed specifically for HPC
- Secure by design
- Native support for GPU acceleration
- **Perfect for medical imaging workloads**
- Supported on our HPC systems!

Basic Container Workflow

1. **Define** your environment
 1. Software versions
 2. Dependencies
 3. Configuration
2. **Build** your container
 1. Create image from definition
 2. Test locally
3. **Run** your analysis
 1. Deploy on HPC
 2. Process your data



Let's think inside
the Box with FSL...

Example: Neuro Image Processing with FSL

```
Bootstrap: docker
From: rockylinux:8.10


%environment
# FSL setup
export FSLDIR=/opt/fsl
export PATH=$FSLDIR/bin:$PATH
export FSLOUTPUTTYPE=NIFTI_GZ
# Python environment
export PYTHONPATH=/opt/scripts:$PYTHONPATH

%post
# System dependencies
dnf -y install epel-release
dnf -y update
dnf -y install \
    python3-pip \
    python3-devel \
    gcc \
    wget \
    mesa-libGL \
    libgomp

# Install FSL
mkdir -p /opt/fsl
wget https://fsl.fmrib.ox.ac.uk/fsldownloads/fslinstaller.py
python3 fslinstaller.py -d /opt/fsl -V 6.0.7.8


# Install Python packages for neuroimaging
python3 -m pip install \
    nibabel \
    nipy \
    scipy \
    numpy
```


Running Containers on HPC

A dark blue terminal window with three colored window control buttons (red, yellow, green) in the top left corner. It contains three lines of white text representing shell commands.

```
$ module load apptainer  
$ cd project/code  
$ apptainer build fsl.sif fsl.def
```

Basic FSL Commands



```
# Run BET on your data
apptainer run \
    --bind $PWD:/data \
    $CONTAINER \
    bet /data/input.nii.gz
    /data/output.nii.gz
```

Key Points to Remember

- **Important Files:**
 - fsl.def: Container definition
 - fsl.sif: Built container
- **Basic Usage**
 - `apptainer run \`
 `--bind /your/data:/data \`
 `/path/to/fsl.sif \`
 `fsl-command`
- **Data Access**
 - Always bind your data directory
 - Use absolute paths
 - Keep data organized

Best Practices for Medical Imaging

1. **Include GPU Support**

1. Essential for deep learning
2. Speeds up image processing

2. **Data Management**

1. Mount data directories
2. Use efficient I/O patterns

3. **Resource Optimization**

1. Request appropriate resources
2. Use parallel processing when possible

Common Use Cases in Medical Imaging

- MRI Analysis Pipelines
- CT Image Processing
- Deep Learning Models
- Image Registration
- Segmentation Workflows

Getting Started

1. **Start Small**
 1. Begin with a simple container
 2. Add complexity gradually
2. **Use Available Resources**
 1. Pre-built containers
 2. HPC documentation
 3. Community support
3. **Document Everything**
 1. Container definitions
 2. Run commands
 3. Workflow steps

Support and Resources

- HPC Help Desk
- Medical Imaging Communities
 - Neuroimaging Tools and Resources (NITRC)
 - BioConda

Thank You!

Remember:

Containers make your research reproducible

Start small and build up

Help is always available

Questions?

Contact: dennist@wustl.edu